

Human ACE2 Protein



Catalog Number: 19CoV-A100

Protein Name: Human ACE2 Protein

Catalog Number: 19CoV-A100

Purity: 95% or above

Alternative Name: Angiotensin-converting enzyme 2, hACE2

Protein Construction: Human ACE2 protein containing Gln18-Asp615 (Accession # [Q9BYF1-1](#)) C-terminal fusion with rabbit IgG Fc.

Endotoxin: <0.1 EU per ug of protein as determined by LAL method

Origin Species: Human, Homo sapiens

Expression Host: HEK293F

Protein type: Recombinant

Purification tag: C-terminal rabbit IgG Fc. Protein A column purified

Molecule Mass: Predicted molecular weight 194 kD

Regulation restrictions: For research use only

Shipping: Recombinant proteins are provided as frozen liquid and will be shipped out with dry ice.

Stability and storage: The protein is stable in a liquid state at -70°C for 12 months. Avoid repeated freeze-thaw cycles.



6440 Lusk Blvd, San Diego, CA 92121

Website: www.exonbio.com

Tel: +1 858-457-1920