

# Recombinant Nucleocapsid, Omicron variant (P13L, E31-, R32-, S33-, R203K, G204R)



**Catalog Number:** 19CoV-S250

---

**Protein Name:** Recombinant Nucleocapsid, Omicron variant (P13L, E31-, R32-, S33-, R203K, G204R)

**Protein Purification:** Full-length SARS-CoV-2 Nucleocapsid protein, Met1-Ala419, with point mutation **P13L, E31-, R32-, S33-, R203K, G204R** and a C-terminal His-tag.

**Alternative Name:** Protein N, NC

**Protein Construction:** N protein with C-terminal His-tag

**Endotoxin:** <0.1 EU/ug

**Expression Host:** E. coli

**Protein type:** Recombinant

**Purification tag:** His-tag

**Formulation:** PBS

**Regulation restrictions:** For research use only

**Shipping:** Recombinant proteins are provided as frozen liquid and will be shipped out with dry ice

**Stability and storage:** The protein is stable in a liquid state at -70°C for 12 months. Avoid repeated freeze-thaw cycles



6440 Lusk Blvd, San Diego, CA 92121

Website: [www.exonbio.com](http://www.exonbio.com)

Tel: +1 858-457-1920