

RSV Nucleoprotein Rabbit Monoclonal Antibody



Product Name: RSV Nucleoprotein Monoclonal Antibody

Immunogen: Influenza B Nucleoprotein

Host: Rabbit

Preparation: Recombinant monoclonal antibody expressed in HEK293F cells and purified by protein A affinity chromatography.

Applications: Any application that detects RSV virus, pseudo virus, and RSV native nucleocapsid protein. Does not detect denatured nucleocapsid protein.

Reactivity: RSV nucleocapsid protein

Isotype: rabbit IgG, kappa Conjugation: Unconjugated Clonality: Monoclonal Concentration: 1mg/ml
Purification: Protein A affinity chromatography

Storage buffer: PBS

Storage conditions: For long term storage, the antibody should be stored at -20°C and is stable for 12 months from the date of receipt. Avoid repeated freeze-thaw cycles. For short term storage, keep at 4°C for a week should be fine.

Application Dilution: 1:1,000-1:100,000 (platform dependent, customers must test the optimal concentration on their own platform)



6440 Lusk Blvd, San Diego, CA 92121

Website: www.exonbio.com

Tel: +1 858-457-1920