

# Rabbit Anti-Ebola Virus VP40 Protein Antibody



---

**Product Name:** Rabbit Anti-Ebola Virus VP40 Protein Antibody

**Immunogen:** VP40

**Host:** Rabbit

**Preparation:** Recombinant monoclonal antibody expressed in HEK293F cells and purified by protein A affinity chromatography.

**Applications:** ELISA, Indirect Immunofluorescence Assay, Lateral Flow

**Reactivity:** Sudan, Mayinga, Bundibugyo and Tai Forest

**Isotype:** rabbit IgG, kappa **Conjugation:** Unconjugated **Clonality:** Monoclonal **Concentration:** 1mg/ml  
**Purification:** Protein A affinity chromatography

**Storage buffer:** PBS

**Storage conditions:** For long term storage, the antibody should be stored at -20°C and is stable for 12 months from the date of receipt. Avoid repeated freeze-thaw cycles. For short term storage, keep at 4°C for a week should be fine.

**Application Dilution:** 1:1,000-1:100,000 (platform dependent, customers must test the optimal concentration on their own platform)



6440 Lusk Blvd, San Diego, CA 92121

Website: [www.exonbio.com](http://www.exonbio.com)

Tel: +1 858-457-1920