

Kappa Free Light Chain (FLC) Rabbit Monoclonal Antibody



Product Name	Kappa Free Light Chain (FLC) rabbit monoclonal antibody
Immunogen	Purified Human Free Light Chain (FLC) - Kappa
Host	Rabbit
Preparation	Recombinant monoclonal antibody expressed in HEK293F cells and purified by protein A affinity chromatography.
Isotype	rabbit IgG, kappa
Conjugation	Unconjugated
Clonality	Monoclonal
Concentration	1mg/ml
Purification	Protein A affinity chromatography
Storage buffer	PBS
Storage conditions	For long term storage, the antibody should be stored at -20°C and is stable for 12 months from the date of receipt. Avoid repeated freeze-thaw cycles. For short term storage, keep at 4°C for a week should be fine.
Application Dilution	1:1,000-1:100,000 (platform dependent, customers must test the optimal concentration on their own platform)



5950 PRIESTLY DR. Ste 101
Carlsbad CA 92008
www.exonbio.com
Tel:(760) 814-8078