

LDH rabbit monoclonal antibody



Product Name	LDH rabbit monoclonal antibody
Immunogen	Plasmodium falciparum (Pf) Plasmodium vivax (Pv) Plasmodium malariae (Pm) Palasmodium ovale (Po)
Host	Rabbit
Preparation	Recombinant monoclonal antibody expressed in HEK293F cells and purified by protein A affinity chromatography.
Isotype	rabbit IgG, kappa
Conjugation	Unconjugated
Clonality	Monoclonal
Concentration	1mg/ml
Purification	Protein A affinity chromatography
Storage buffer	PBS
Storage conditions	For long term storage, the antibody should be stored at -20°C and is stable for 12 months from the date of receipt. Avoid repeated freeze-thaw cycles. For short term storage, keep at 4oC for a week should be fine.
Application Dilution	1:1,000-1:100,000 (platform dependent, customers must test the optimal concentration on their own platform)

